



# ESF5-18 SERIE

## SURGICAL ASPIRATORS



Portable suction machine for minor surgery and post-tracheotomy, for home use.  
Ideal for home post-surgery treatment.  
Anti-shock case, light, easy to carry and keep clean.  
Supplied with 1 lt. jar sterilizable, with over flow valve.  
Complete with suction regulator, vacuum gauge.  
Filter disposable

### TECHNICAL SPECIFICATIONS

- Suction power : 20 lt./min.
- Maximum vacuum : 750 mm/hg.
- Power consumption: 90 W
- One safety fuse
- Voltage: 230V - 50/60 Hz
- Size: 260 x 180 x 220 (H) mm
- Weight: 3,4 kg

Available also with jar of 2 lt capacity

Aspirateur chirurgie mineure et post-trachéotomie à utiliser a domicile.  
Idéal pour traitements post-opératoires à domicile.  
Boîtier en ABS antichoc, léger et facile à nettoyer.  
Livré avec bocal stérilisable de 1 l., avec soupape de trop-plein.  
Régulateur d'aspiration avec vacuomètre.  
Filtre usage unique,

### CARACTÉRISTIQUES TECHNIQUES

- Puissance d'aspiration: 20 lt./min.
- Aspiration maximum: 750 mm/hg.
- Puissance absorbée: 90 W
- Un fusible de sûreté.
- Alimentation: 230V - 50/60 Hz
- Dimensions: 260 x 180 x 220 (H) mm
- Poids: 3,4kg

Disponibile aussi avec bouteille 2 lt



ESSE 3 Via Garibaldi 30  
14022Castelnuovo D.B. (AT)  
tel +39 011 99 27 706  
fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web : [www.esse3-medical.com](http://www.esse3-medical.com)

