




ERS-1900P

Flat screen chart monitor



OPHTHALMOLOGY / ENT




ERS-1900P

Flat Screen Chart Monitor



- Polarized chart for binocular test
- 23-inch high contrast and bright display
- Smooth and quiet operation
- Wide range of test types
- Precise distance adjustment
- Easy mounting
- Random test sequence



Rotating Astigmatic Fan Chart

Fast accurate test for cylinder and axis



Computerized vision tester



Support data connection of the ERS-1900P with a computerized vision tester



Technical Features

Display	23" LCD (color) High resolution Monitor
Brightness	250 cd/m ²
Chart versions	E, C, ABC, 123, child, special version charts
Chart acuity	6/3-6/150 20/10-20/500 0.04-2.0
Distance	2.0 ~ 7.0m
Power supply	Input: AC100~240V ,50/60Hz
Power consumption	<=33W
Dimension (WxHxD)	530x390x120 mm
Weight	5.5kg (Net. Weight)
Remote	IR
Batteries for remote control	3A



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

