



DATA SHEET 

TRIAL FRAMES

Cod. **RS-400** Trial Frame

PD Scale: 50mm~80mm
Net weight: 61g
Knobs fastened with screws
Back lens holders of round corners
Easily-shaped ear rests
Available for 5 pieces trial lensed to be inserted at each side



Cod. **RS-300** Trial Frame

Titanium alloy trial frame
P.D. is adjustable from 54mm to 70mm.
Axis of astigmatism.
One touch adjustment on nose pad
Easy slide-on temples
Temple with spring hinge.
Available in angle and length
Accommodate up to four 38mm trial lenses.
light 47g, simple ,reasonable



Cod. **RS-100** Trial Frame

Plastics , PD adjustable. PD 48-80mm

....



Cod. **RS-200** Trial Frame

Plastics frame. PD is not adjustable.
packaged in one box, one box with five pcs of frames.
PD: 48# 50# 52# 56#,58#,60#62#64#66# 68# 70#



Cod. RS-500A Trial Frame (Child)

Plastic material ,
for child,
PD size: 48mm, 56mm, two type



Cod. RS-600A Trial Frame

New style, high quality. There are PD
48 50 54 56 58 60 62 64 66 68 70mm.



Trial Clip

With or without gradient



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

