



# STREAK RETINOSCOPE

**Rechargeable**



**Cod. YZ24B**



**Optional:**



- The filament can be 360°rotated and upward and downward moved.
- Quickly and precisely measure the astigmatism axis.
- The light can be converged, radiated and paralleled.
- Continuously adjustable brightness.
- Convenient quick recharging.
- Anti-recharging and anti-discharging function protects the battery.
- Can also be connected with the head of ophthalmoscope.

## TECHNICAL FEATURES

Cod.	<b>YZ24B</b>
Working Distance	0.67m
Streak	3mm-20mm
Illumination Source	3.5V/2.8W halogen bulb
Streak Rotation	360°
Input Voltage	AC220V±22V/50Hz±1Hz, AC110V±11V/60Hz±1Hz
Power(Main Body of Ophthalmoscope)	Li+ battery pack DC:3.7V



ESSE3 srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3-medical.com](http://www.esse3-medical.com)

