




APS-T00

Projection perimeter



OPHTHALMOLOGY / ENT



APS!HCC

..... Dfc^YVW|cb`DYf]a YHf




Cod. 5DG!HCC

- Full-field projection perimeter
- Software and hardware designed in line with international standards.
- Automated, user-friendly design, easier operation
- Automatic pupil measurement, monitoring and correction of eye position
- 3 minutes to complete the threshold perimetry
- Report analysis is consistent with the international mainstream.

Technical Features

Device Type	Integrated Type
Stimulate point appearing way	Optical projection
Detection method	Static perimetry
Detection range	0-90 degree
Background light	Yellow , white
Background brightness	White:31.5asb Yellow:315asb
Visual Target color	Red,blue,white
Visual Target Brightness	Oasb-10000asb
Vision Target intervaltime	≥200ms(customizable,no limit)
Visual target holding time	≥200ms(customizable,no limit)
Brightness level	0-51db
Visual Target Size(Goldman)	I Grade II Grade III Grade V Grade
Detection strategy	Full-threshold, fast threshold, custom, Two notation, three notation, quantify defects
Threshold Policy	Intelligent dynamic, fast intelligent dynamic, age-related, The threshold, a single strength
Threshold Test Mode	Center 10-2, Center 24-2, Center 30-2, Center 60-4, Nasal step, macula
Fixation monitoring	Dynamic real-time video monitoring, physiological blind spot monitoring, Eye tracking curve, eye position offset alarm
Analysis Software	Reliability analysis, analysis of a single vision, triple report analyzes
Report Figure	Value, grayscale, the overall bias decibel chart, The overall deviation probability map, chart pattern deviation decibels, Pattern deviation probability plot, staring Figure



Full Field Dynamic test Report

Name: DynamicDemo Sex: Male Age: 29 Eye: R ExamID: 20141022000


A:20.4 deg III White

Pupil Diameter:0

Date:2014-10-22

Time:10:51:18

Age:29



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

