



# ES-FA6500 ES-FA6500K

Auto Refractometer



OPHTHALMOLOGY / ENT



## Features

ES-FA6500K with keratometer function.  
ES-FA6500 Without keratometer function, others same.

### Technical Data

- Touch screen, more easy operate.
- Tilt screen.
- Motorized chinrest.
- With keratometer, auto tracking function.



Sphere	-20~+20D(VD=12)0.125D step
Cylinder	-8~+8D,0.125D Step
Axis	0~180°1°step
PD	45~88mm, 1mm step
VD	0mm, 12mm, 13.75mm
Min.Pupil Size	2mm
Radius of curvature	5.0~10mm(increment:0.01mm)
Corneal Power	33.75-67.50D (when corner equivalent refractive index is 1.337) (increment selectable from 0.12, 0.25D)
Corneal astigmatism	0~8D(increment selectable from 0.12, 0.25D)
Axis	1~180°(increment:1°)
Corneal diameter	2.0~14 mm(increment:0.1mm)
Memory of data	10 measure vature for each right and left eye
Chart	Follow up colorful fogging chart system
display	7" TFT LCD
print	Thermal printer
Power supply	220V/110V, 50/60Hz
Dimension	478*268*472mm
Weight	-20kg(44lb)



ESSE3 srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3.dreamgest.com](http://www.esse3.dreamgest.com)

