



ES-HE100

Portable Keratometer



OPHTHALMOLOGY / ENT



ES-HE100

Portable Keratometer

Technical specification

Keratometer ES-HE100: Electronic and optical integration, ultra small corneal curvature tester, is mainly used to measure the corneal curvature radius and diopter, wireless output print data.

Measuring range: 6.5mm~9.5mm

Precision: +/-0.05mm

Resolution of curvature radius of cornea: 0.01mm

single measuring time: 0.03s

Output: wireless infrared thermal printer

Can observe the eye directly through the screen.

Weight: <0.5Kg(with batteries)

Dimension: 240mmX90mmX60mm

Power: 500mW+15%



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

