



# JD-2600A & JD-2600B

## Auto-Lensmeters



OPHTHALMOLOGY / ENT



# JB-2006A & JB2006B

## Auto-Lensmeters



- Cod. **JD-2006A** - Auto Lensmeter with 7" screen
- Cod. **JD-2006B** - Auto Lensmeter with 5,6" screen

### Adjustable LCD screen

Tilt color screen provides you with clear view from any position, such as sitting and standing.



### Sensitive touch LCD display screen with pen



Just touch the screen can finish your measurement, your operation could be easier and more convenient. This new generation auto lens meter is based on Hartmann Sensor technology. Advanced simultaneous measurement of 108 points provides a faster measurement with greater accuracy and reliability.

### Measurement screen

Easy to use with recognisable icons. The JD-2600A and JD-2600B provide a user-friendly graphical interface for an immediate understanding.



### Automatic recognition of single and progressive measurement



Automatic measurement of progressive lenses. A graph helps the operator to improve accuracy and speed in the measurement of near and far vision.

### Built-in UV tester and Hi-speed UV measurement

Just hold the lens through the sensor, the display will easily indicate the UV measurement by the progress bar



### Built-in thermal printer



Quick data print with built-in printer. Easy to reload printing paper.



## Technical Features

Range of measurement	
<b>Sphere</b>	±25D 0.01/0.12/0.25D step
<b>Cylinder</b>	±9.99D 0.01/0.12/0.25D step
<b>Axis</b>	0°~180°(1° step)
<b>Add.</b>	0~+9.99D 0.01/0.12/0.25D step
<b>Prism degree</b>	0~15Δ 0.01/0.12/0.25Δ step

Measurement mode	
<b>Cylinder</b>	+,+/-,-
<b>Prism</b>	X-Y,P-B
<b>Contact lens</b>	Soft/hard
<b>Measuring mode</b>	Single/progressive/automatic recognition

Specifications	
<b>Lens diameter</b>	20 ~ 108mm
<b>PD</b>	40~90mm, 0.5mm step
<b>Measurement speed</b>	0.1s
<b>Display</b>	7" TFT LCD (Cod. <b>JD-2006A</b> ) 5.6" TFT LCD (Cod. <b>JD-2006B</b> )
<b>Printer</b>	Thermal printer
<b>Dimension</b>	220x320x430 mm
<b>Weight</b>	Approx. 7KG
<b>Power supply</b>	100~240V 50~60HZ



ESSE3 srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3-medical.com](http://www.esse3-medical.com)

