



DATA SHEET



Model: ESR-B

OPHTHALMIC/ENT CHAIR



Description

With tiltable backrest.
Can move up and down by motor.
Maintenance free linear actuator lift system
Chair elevation control with the footswitch
Durable, stainless steel kick plate
Upholstered arms that rotate upward
Durable, easy-to-clean upholstery
Height to headrest: 128 cm
Depth when fully tilted: 128 cm
Seat Dimensions: 53 cm x 48 cm
Chair Elevation: 57 cm to 78 cm



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

