

BHE-A



Product Data **Height adjustable bucky table**

HORIZONTAL BUCKY TABLE WITH ELEVATING TABLE TOP

MAIN FEATURES

- Four-way floating patient table top with motorized height adjustment.
Longitudinal travel of the under-table Potter Bucky carriage with electromagnetic brakes control.
Foot pedals control for up/down motorized motion and table top electromagnetic brakes.

TECHNICAL DATA

HORIZONTAL TABLE

- Table top dimensions 2200 x 810 mm
- Longitudinal travel 900 mm (\pm 450 mm)
- Transversal travel 220 mm (\pm 110 mm)
- Table top - floor distance Minimum 500 mm / maximum 825 mm
- Table top vertical travel Motorized
- Maximum patient's weight 200 Kg
- Brakes Electromagnetic, foot-controlled
- Weight 220 Kg
- Power supply 230Vac (\pm 10%) single phase - 50/60 Hz

POTTER BUCKY

- For all standard radiographic cassettes Size from 13x18 cm up to 35x43 cm, both direction
- Longitudinal travel 550 mm (\pm 275 mm)
- Grid Vibrating R 10:1 F 120 cm 103 lines /inch
- Brakes Electromagnetic, push button-controlled

OPTIONALS

- Compression band
- Handgrips
- Lateral cassette holder




BHE-A

Product Data **Height adjustable bucky table**



OVERALL DIMENSIONS (mm)



PD-127-E 02/11/2015 Rev.00

Data subject to modification without prior notice



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3.dreamgest.com

