



ES-DTR Series

CHEST BUCKY



IMAGING SYSTEMS



ES-DTR Series

Commercial Name	ES-DTR - Chest Bucky	
REF	ES-DTRXXX	
Class	The medical device is in class I	
Product family	CODE	DESCRIPTION
	ES-DTR600BVS	Chest stand designed for wall fixation, panel in ABS with filtration 0.5mm Al.Eq.. Complete of Potter Bucky for radiographic cassette size up to 430x430mm, set for AEC. Oscillating anti-scattering grid size 480x438mm R10:1 F150 103lp/inch not included. Accessories for patient positioning available at your choosing.
	ES-DTR604BVS	Chest stand designed for wall fixation, panel in ABS with filtration 0.5mm Al.Eq.. Complete of digital detector support for cassette size up to 430x430x17mm (17"x17"x0,6") with support for removable anti-scattering grid size 480x438mm, set for AEC. Grid not included. Accessories for patient positioning available at your choosing.
	ES-DTR604500BVS	Chest stand designed for wall fixation or for fixation with self supporting base included, panel in ABS with filtration 0.5mm Al.Eq.. Complete of digital detector support for cassette size up to 430x430x17mm (17"x17"x0,6") with support for removable anti-scattering grid size 480x438mm, set for AEC. Grid not included. Accessories for patient positioning available at your choosing.

TECHNICAL FEATURES

ES-DTR600BVS

Features and overall dimensions:


- Front panel: ABS
- Filtration: 0.5 Al eq.
- Stand fixation: on wall
- Vertical Potter-Bucky movement: 1565 mm
- Minimum distance from floor of imaging area: 385 mm
- Maximum distance from floor of imaging area: 1950 mm
- Vertical movement balance: Counterweight
- Vertical movement brake: mechanical
- Power supply: 24Vdc 11VA
- Standard Potter-Bucky grid: R 10 103 lines per inches
S.I.D. 150 cm
- Device Weight: 105 kg

ES-DTR604BVS/ES-DTR604500BVS

Features and overall dimensions:


- Front panel: ABS
- Filtration: 0.5 Al eq.
- Stand fixation: on wall
- Vertical detector movement: 1565 mm
- Minimum distance from floor of imaging area: 385 MM
- Maximum distance from floor of imaging area: 1950 mm
- Vertical movement balance: Counterweight
- Vertical movement brake: mechanical
- Device Weight: 90 kg

ES-DTR Series





ES-DTR600BVS

All the following measures are in mm



ES-DTR604BVS/ES-DTR604500BVS

All the following measures are in mm



Accessories

- cod. AS045 – Self-support base
- cod. AS030 – Arm support
- cod. AS020 – Lateral handles
- cod. CL0001 - Varex AEC solid state measurement chamber 3fields with preamplifier and 20 mt extension cables for generator with WECO connector.
- cod. CL0002 - Varex AEC solid state measurement chamber 3fields with 20 mt extension cables for generator with WECO connector (no preamplifier).
- cod. CL0003 - Varex AEC solid state measurement chamber 3fields with preamplifier and 15 mt extension cables for generator with LEMO connector.
- cod. CL0004 - Varex AEC solid state measurement chamber 3fields with 15 mt extension cables for generator with LEMO connector (no preamplifier).
- cod. M0046 – Grid 215L” R:10 S.I.D.150cm 480x438mm.

