




DATA SHEET



Model: **BBR 925 xPRO**

Blood Bank Refrigerator



Description

Color: White

Panel: xPRO (BBR)

- 3.5" TFT color control panel with capacitive buttons and password
- Automatic test facilities for temperature alarms and functioning status
- Backup battery for managing alarms in case of power failure
- Audible and visual alarms
- Temperature chart recorder

Technical Features

- Temperature +4°C
- Refrigeration Type Forced air
- Defrost Type Automatic, with heating resistance
- Condensate Evaporation Automatic
- Refrigerant Gas R290
- Climatic Class +16°C / +38°C (ST)
- External Dimensions (W x D x H mm) 1200 x 770 x 1950

Electrical Data

- Power Supply (V/ph/Hz) 220-230/1/50
- Power Plug Type F

Structural Features

- Structure: Upright
- External Material: Anti-bacterial coated steel white colour
- Internal Material: Stainless steel 18/10 AISI 304
- Insulation Thickness (mm): 50
- Triple-Layer Glass Door :Yes
- Automatic Door Closure: Yes
- Drawers with extractable slides: 12
- Drawer Material: Stainless steel 18/10 AISI 304
- Plastic Cross Dividers: 6 each drawer
- Blood Bag Storage Capacity (450 ml): 360

Control Panel


- Display Type: TFT 3,5"
- Languages: ITA - ENG - FRA
- Access Password: Yes
- Temperature Graph Visible on display
- Probe Type: Nr. 3 PT1000
- Alarm Signaling: Acoustic and visual
- Alarm Indications: Text form / codes
- Adjustable High/Low Temperature Alarm: Yes
- Adjustable Door Open Alarm: Yes
- Power Loss Alarm: Yes
- Probe Fault Alarm: Yes
- Low Battery Alarm: Yes
- Dirty Condenser Alarm: Yes
- Evaporator Freezing Alarm: Yes
- Acoustic Alarm Silencing Options: Yes
- Acoustic Alarm Reactivation: Yes

- Alarm Functionality Test: Yes
- Calibration: Yes
- Alarm Memory: 30
- Door Openings Memory: 30 days
- Defrost Memory: 30 days
- USB Port: Optional
- Data and Event Storage: Yes (on internal memory)
- Backup Battery for Display and Alarms: Yes
- Backup Battery Autonomy for Display and Alarms (hours) : 48

ACCESSORIES AVAILABLE



Plastic partitions for drawers



Front top for drawers



Kit of castors




Safety door lock with key




Data logger function with USB port




Ballasted sensor for simulation of product temperature alarm



Cable port with cover Ø 18 mm



Cable port with cover Ø 45 mm



Back-up cooling unit




Deduction for solid door (instead of glass door)



Air sterilization lamp



Safety anti-freezing thermostat



Temperature uniformity and stability kit



Natural refrigerant R290




GSM/GPRS telephone dialer



Different voltage



Different voltage



Wooden crate packing



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3.dreamgest.com

