




ESPS-100AD

Digital Ultrasonic Cleaner
Double Frequency Series



STERILIZATION



ESPS-100AD

Digital Ultrasonic Cleaner Double Frequency Series

SPECIFICATIONS

Model ESPS-100AD
Power 600W
Capacity 30L
Frequency 28/40KHz
Tank size 500*300*200mm (L×W×H)
Overall size 550*330*360 mm (L×W×H)
Set time 1-99minutes
Set temperature normal temperature-80
Heating power 500W
Power supply AC220V 50Hz

FEATURES

Service line, Electron line, Medicament line, Laboratory, Machine line, Horniness alloy line, Chemistry line, etc

ACCESSORISE:

- Machine
- Cover
- Power line
- Name Plate
- Certificate



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

